

SpiroGuide

Heads-Up Display

- ① Hands-free, at-a-glance cylinder pressure information
→ **Allows concentration on the situation at hand**
- ② Highly visible Buddy Reading
→ **All team members can stay informed about each other's status**
- ③ Installation inside the mask
→ **Optimal visibility and protection**
- ④ Wireless technology
→ **No cables. Minimal installation time.**
- ⑤ Two-way wireless communication
→ **100% secure link and no risk for mismatch between HUD and BAC**
- ⑥ Turn-Back signal indication
→ **Get out in time**
- ⑦ Installation without the need for tools
→ **Easy to install and remove**
- ⑧ Status indication of ADSU / PASS
→ **Warnings are both auditory and visual**



SpiroGuide

Heads-Up Display

The lights in the **Heads-Up Display (HUD)** instantly and effectively show the cylinder pressure of your equipment while also showing your colleagues your status. Whether it is in heavy smoke or bright daylight – the lights adapt to the conditions as you concentrate on the situation at hand. The light indications can be tailor-made to suit all operational procedures.



Thanks to SpiroGuide's wireless two-way communication, you get a 100% secure link between the **Breathing Apparatus Computer (BAC)** and the wireless Heads-Up Display. Furthermore you never risk a mismatch. You can be sure your HUD will always display the information from your BAC. And going wireless means no cables and minimal installation time.

An optional **Turn-Back signal** ensures you get out in time. The level can either be set automatically or to suit your own operating procedures.

A flashing LED on the right-hand side informs you that your unit has entered **ADSU / PASS** pre-alarm. If the LED is permanently lit, the unit is in full alarm mode. In this way, you stay informed about the ADSU / PASS status at all times, even in very noisy environments.

The Heads-Up Display is equipped with a **"Buddy Reading"** feature. Your colleagues can easily read the lights on your HUD, just as you can read theirs. The SpiroGuide Heads-Up Display is the reliable way for you to stay in control at all times. It's also the best solution for the team.



STANDARDS

- EN50020
- EN50014
- EN50082-1
- EN61000-6
- UL913
- JCDD/38
- NFPA 1981
- EN137:2006

APPROVALS

- Directive PPE 89/686/EEC
- Directive EMC 89/336/EEC
- NIOSH 42 CFR (Part 84)
- NFPA 1981:2002

BATTERY DATA

- 2 x AAA alkaline batteries
- Duracell MN2400 or Energizer E92 recommended for optimum performance
- Up to 50 hours during normal use

WEIGHT

- 55 gr including batteries

OPTIONS

- Turn-Back signal
- ADSU / PASS status indication

INTERSPIRO

www.interspiro.com

CENTRAL EUROPE

AUSTRIA

INTERSPIRO GesmbH www.interspiro.de
Sternweg 20 A-8401 ZETTLING AUSTRIA
TEL +43 (0)313 557 333 FAX +43 (0)313 557 333 22 E-MAIL info@interspiro.de

GERMANY

INTERSPIRO GmbH www.interspiro.de
Postfach 1220 D-76691 FORST/BADEN GERMANY
TEL +49 (0)7251 8030 FAX +49 (0)7251 2298 E-MAIL info@interspiro.de

SWITZERLAND

INTERSPIRO AG www.interspiro.de
Güterstraße 47 CH-4133 PRATTELN SWITZERLAND
TEL +41 61 827 99 77 FAX +41 61 827 99 70 E-MAIL info@interspiro.ch

THE NETHERLANDS & BELGIUM

INTERSPIRO BV www.interspiro.nl
Televisieweg 113 NL-1322 BD ALMERE NETHERLANDS
TEL +31 (0)36 5363103 FAX +31 (0)36 5384809 E-MAIL infobv@interspiro.com

NORTH & SOUTH AMERICA

INTERSPIRO Inc. www.interspiro-us.com
10225 82nd Avenue PLEASANT PRAIRIE WI 53158-5801 USA
TEL +1 262 947 9000 FAX +1 262 947 9020 E-MAIL info@interspiro-us.com

UNITED KINGDOM & IRELAND

INTERSPIRO Ltd. www.interspiro.com
7 Hawksworth Road Central Park TELFORD Shropshire TF2 9TU UNITED KINGDOM
TEL +44 (0)1952 200 190 FAX +44 (0)1952 299 805 E-MAIL info@interspiro.com

SCANDINAVIA, ASIA/PACIFIC & MIDDLE EAST

SWEDEN

NORDIC & EXPORT SALES DIVISION www.interspiro.com
Box 10060 S-181 10 LIDINGÖ SWEDEN
TEL +46 8 636 51 00 FAX +46 8 765 48 53 E-MAIL info@interspiro.com

MALAYSIA

NORDIC & EXPORT SALES DIVISION www.interspiro.com
305 & 305A Lorong Perak Taman Melawati 53100 KUALA LUMPUR MALAYSIA
TEL +60 3 4105 8122 FAX +60 3 4105 3122 E-MAIL asiapacific@interspiro.com